



Product Datasheet

EAF1 antibody (orb394995)

Description	Rabbit polyclonal antibody to EAF1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	ELL associated factor 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	43 kDa
Uniprot ID	Q96JC9
Dilution Range	WB: 1:500-1:1000; IP: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K-562 cells using EAF1 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using EAF1 antibody.