



Product Datasheet

ECHS1 antibody (orb394996)

Description	Mouse monoclonal antibody to ECHS1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	enoyl Coenzyme A hydratase, short chain, 1, mitochondrial
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	31 kDa
Uniprot ID	P30084
Dilution Range	WB: 1:2000-1:20000; IP: 1:500-1:2000; IHC: 1:50-1:1000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using ECHS1 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using ECHS1 antibody.