



Product Datasheet

EEF1D antibody (orb395000)

Description	Mouse monoclonal antibody to EEF1D.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	eukaryotic translation elongation factor 1 delta(guanine nucleotide exchange protein)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	40 kDa
Uniprot ID	P29692
Dilution Range	WB: 1:500-1:2000; IF: 1:10-1:100; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical staining of paraffin...