


Product Datasheet

EIF2C3-Specific antibody (orb395005)

Description	Rabbit polyclonal antibody to EIF2C3-Specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	eukaryotic translation initiation factor 2C, 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	90-100 kDa
Uniprot ID	Q9H9G7
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of SKOV-3 cells US...