



Product Datasheet

EIF3D antibody (orb395006)

Description	Mouse monoclonal antibody to EIF3D.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	eukaryotic translation initiation factor 3, subunit D
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	66 kDa
Uniprot ID	O15371
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using EIF3D antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using EIF3D antibody.