


Product Datasheet

EIF3M antibody (orb395007)

Description	Mouse monoclonal antibody to EIF3M.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	eukaryotic translation initiation factor 3, subunit M
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	43 kDa
Uniprot ID	Q7L2H7
Dilution Range	WB: 1:500-1:2000; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using...