



Product Datasheet

ESR2 Beta-6 antibody (orb395017)

Description	Rabbit polyclonal antibody to ESR2 Beta-6.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	estrogen receptor 2(ER beta)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	42 kDa
Uniprot ID	Q92731
Dilution Range	WB 1:500-1:1000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse spleen ti...



Immunohistochemical staining of paraffin...