



Product Datasheet

Factor XII antibody (orb395022)

Description	Mouse monoclonal antibody to Factor XII.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	coagulation factor XII(Hageman factor)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	28 kDa
Uniprot ID	P00748
Dilution Range	WB: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of K-562 cells using Factor XII antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using Factor XII antibody.