


Product Datasheet

FADS antibody (orb395023)

Description	Rabbit polyclonal antibody to FADS.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	FAD1 flavin adenine dinucleotide synthetase homolog(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa, 65-70 kDa
Uniprot ID	Q8NFF5
Dilution Range	WB: 1:500-1:1000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K-562 cells using FADS antibody (orb395023).