



Product Datasheet

FAM175A antibody (orb395026)

Description	Rabbit polyclonal antibody to FAM175A.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	family with sequence similarity 175, member A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	46 kDa
Uniprot ID	Q6UWZ7
Dilution Range	WB: 1:500-1:2400; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human testis ti...



Immunohistochemical staining of paraffin...