


Product Datasheet

FAM19A5 antibody (orb395028)

Description	Rabbit polyclonal antibody to FAM19A5.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	family with sequence similarity 19(chemokine(C-C motif)-like), member A5
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	14 kDa
Uniprot ID	Q7Z5A7
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...