



Product Datasheet

FAM8A1 antibody (orb395033)

Description	Rabbit polyclonal antibody to FAM8A1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	family with sequence similarity 8, member A1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	44 kDa
Uniprot ID	Q9UBU6
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK-293 cells u...



Immunohistochemical staining of paraffin...