



Product Datasheet

FANCD2 (Phospho-S330) antibody (orb395035)

Description	Rabbit polyclonal antibody to FANCD2 (Phospho-S330).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fanconi anemia, complementation group D2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	150-155 kDa
Uniprot ID	Q9BXW9
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HL-60 cells using FANCD2 (Phospho-S330) antibody (orb395035).



Immunohistochemical staining of paraffin-embedded tissue sections using FANCD2 (Phospho-S330) antibody (orb395035).