



Product Datasheet

Fc epsilon RI alpha antibody (orb395037)

Description	Rabbit polyclonal antibody to Fc epsilon RI alpha.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Fc fragment of IgE, high affinity I, receptor for; alpha polypeptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	55 kDa
Uniprot ID	P12319
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of human spleen ti...



Immunohistochemical staining of paraffin...