


Product Datasheet

Fetuin-A antibody (orb395040)

Description	Mouse monoclonal antibody to Fetuin-A.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	alpha-2-HS-glycoprotein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	58 kDa
Uniprot ID	P02765
Dilution Range	WB: 1:2000-1:10000; IP: 1:1000-1:6000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using Fetuin-A antibody (orb395040).