



Product Datasheet

fibrillin 2-Specific antibody (orb395043)

Description	Rabbit polyclonal antibody to fibrillin 2-Specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	fibrillin 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	290 kDa, 160 kDa
Uniprot ID	P35556
Dilution Range	WB: 1:500-1:2000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical staining of paraffin...