



Product Datasheet

FKBP2 antibody (orb395047)

Description	Mouse monoclonal antibody to FKBP2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	FK506 binding protein 2, 13kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	16kd
Uniprot ID	P26885
Dilution Range	WB: 1:500-1:1000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human cerebellu...



Immunohistochemical staining of paraffin...