



Product Datasheet

TLS antibody (orb395058)

Description	Mouse monoclonal antibody to TLS.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	fusion(involved in t(12;16) in malignant liposarcoma)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	73 kDa
Uniprot ID	P35637
Dilution Range	WB: 1:500-1:2000; IP: 1:5000-1:50000; IHC: 1:500-1:2500; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of K-562 cells using TLS antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using TLS antibody.