



Product Datasheet

GDI2 antibody (orb395075)

Description	Mouse monoclonal antibody to GDI2.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	GDP dissociation inhibitor 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	46-50kd
Uniprot ID	P50395
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using GDI2 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using GDI2 antibody.