



Product Datasheet

GEFT-Specific antibody (orb395077)

Description	Rabbit polyclonal antibody to GEFT-Specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	RhoA/RAC/CDC42 exchange factor
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	64 kDa
Uniprot ID	Q86VW2
Dilution Range	WB: 1:500-1:3000; IP: 1:500-1:3000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse skeletal ...



Immunohistochemical staining of paraffin...