



Product Datasheet

GLAST antibody (orb395082)

Description	Rabbit polyclonal antibody to GLAST.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Excitatory amino acid transporter 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50-55, 90-100 kDa
Uniprot ID	P43003
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Neuro-2a cells ...



Immunohistochemical staining of paraffin...