



Product Datasheet

GLI3-Specific antibody (orb395084)

Description	Rabbit polyclonal antibody to GLI3-Specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	GLI family zinc finger 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	230kd, 83-86kd
Uniprot ID	P10071
Dilution Range	WB: 1:200-1:1000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of human placenta ...



Immunohistochemical staining of paraffin...