


Product Datasheet

GSTT2B antibody (orb395105)

Description	Rabbit polyclonal antibody to GSTT2B.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	glutathione S-transferase theta 2B(gene/pseudogene)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	28kd
Uniprot ID	P0CG30
Dilution Range	WB: 1:500-1:3000; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using GSTT2B antibody (orb395105).