




Product Datasheet

IL1RAPL2 antibody (orb39511)

Description	Rabbit polyclonal antibody to IL1RAPL2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	IL1RAPL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78670
Uniprot ID	Q9NP60
NCBI	NP_059112.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in K562 cell li...



Western blot analysis of in mouse NIH-3T...