

Product Datasheet

HLA class I (HLA-A) antibody (orb395139)

Description Rabbit polyclonal antibody to HLA class I (HLA-A).

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen major histocompatibility complex, class I, A

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 44 kDa


Uniprot ID **P30443**, **P04439**

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200

Expiration Date 12 months from date of receipt.



Western blot analysis of HepG2 cells using HLA class I (HLA-A) antibody (orb395139).



Immunohistochemical staining of paraffin-embedded tissue sections using HLA class I (HLA-A) antibody (orb395139).