


Product Datasheet

HMGXB4 antibody (orb395145)

Description	Rabbit polyclonal antibody to HMGXB4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	HMG box domain containing 4
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	66 kDa
Uniprot ID	Q9UGU5
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using...