


Product Datasheet

IFITM2-specific antibody (orb395166)

Description	Mouse monoclonal antibody to IFITM2-specific.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	interferon induced transmembrane protein 2(1-8D)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	15 kDa
Uniprot ID	Q01629
Dilution Range	WB: 1:500 - 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using the IFITM2-specific antibody (orb395166).