


Product Datasheet

IFT25 antibody (orb395168)

Description	Rabbit polyclonal antibody to IFT25.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	heat shock protein family B(small), member 11
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	25 kDa
Uniprot ID	Q9Y547
Dilution Range	WB: 1:500 - 1:2400;IF: 1:20 - 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using IFT25 antibody.