



Product Datasheet

IL1F8 antibody (orb395183)

Description	Mouse monoclonal antibody to IL1F8.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	interleukin 1 family, member 8(eta)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	44 kDa
Uniprot ID	Q9NZH7
Dilution Range	WB: 1:2000-1:16000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Recombinant pro...



Immunohistochemical staining of paraffin...