


Product Datasheet

ING1-specific antibody (orb395193)

Description	Rabbit polyclonal antibody to ING1-specific.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	inhibitor of growth family, member 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa
Uniprot ID	Q9UK53
Dilution Range	WB: 1:500-1:1000 IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...