



Product Datasheet

iPLA2 antibody (orb395211)

Description	Rabbit polyclonal antibody to iPLA2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	phospholipase A2, group VI(cytosolic, calcium-independent)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	80 kDa, 88 kDa
Uniprot ID	O60733
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse testis ti...



Immunohistochemical staining of paraffin...