



Product Datasheet

JP45 antibody (orb395217)

Description	Rabbit polyclonal antibody to JP45.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	junctional sarcoplasmic reticulum protein 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30kd
Uniprot ID	Q96MG2
Dilution Range	WB: 1:500-1:1000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human skeletal ...



Immunohistochemical staining of paraffin...