



Product Datasheet

KDM3A,JMJD1A antibody (orb395225)

Description	Rabbit polyclonal antibody to KDM3A,JMJD1A.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	lysine(K)-specific demethylase 3A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	150 kDa
Uniprot ID	Q9Y4C1
Dilution Range	WB: 1:1000-1:4000; IP: 1:500-1:2000 ; IHC: 1:20-1:200; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using JMJD1A antibody.



JMJD1A antibody