



Product Datasheet

Kindlin 2 antibody (orb395241)

Description	Rabbit polyclonal antibody to Kindlin 2.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	fermitin family homolog 2(Drosophila)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	78 kDa
Uniprot ID	Q96AC1
Dilution Range	WB: 1:500-1:3000; IP: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 cells using Kindlin 2 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using Kindlin 2 antibody.