



Product Datasheet

LC3B-Specific antibody (orb395265)

Description	Rabbit polyclonal antibody to LC3B-Specific.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	microtubule-associated protein 1 light chain 3 beta
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	15 kDa
Uniprot ID	Q9GZQ8
Dilution Range	WB:1:500-1:5000, IHC:1:20-1:200, IF/ICC:1:20-1:200, FACS:N/A
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain tissue...



Immunohistochemical staining of paraffin...