



Product Datasheet

LGALS3BP antibody (orb395269)

Description	Mouse monoclonal antibody to LGALS3BP.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	lectin, galactoside-binding, soluble, 3 binding protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	65-70 kDa
Uniprot ID	Q08380
Dilution Range	WB: 1:1000-1:4000; IP: 1:500-1:2000; IHC: 1:50-1:500; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human blood tis...



Immunohistochemical staining of paraffin...