



Product Datasheet

Loricrin antibody (orb395280)

Description	Rabbit polyclonal antibody to Loricrin.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	loricrin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	26 kDa
Uniprot ID	P23490
Dilution Range	WB: 1:500-1:10000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cells using Loricrin antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using Loricrin antibody.