


Product Datasheet

L-Plastin antibody (orb395283)

Description	Rabbit polyclonal antibody to L-Plastin.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	lymphocyte cytosolic protein 1(L-plastin)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	64-70 kDa
Uniprot ID	P13796
Dilution Range	WB: 1:1000-1:6000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cells usin...