


Product Datasheet

MKS3 antibody (orb395311)

Description	Rabbit polyclonal antibody to MKS3.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	transmembrane protein 67
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q5HYA8
Dilution Range	IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...