



Product Datasheet

Myosin Light Chain 2 antibody (orb395332)

Description	Mouse monoclonal antibody to Myosin Light Chain 2.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	myosin, light chain 2, regulatory, cardiac, slow
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	19 kDa
Uniprot ID	P10916
Dilution Range	WB: 1:1000-1:5000; IHC: 1:100-1:1000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart tis...



Immunohistochemical staining of paraffin...