



Product Datasheet

NDUFS5 antibody (orb395342)

Description	Mouse monoclonal antibody to NDUFS5.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	NADH dehydrogenase(ubiquinone) Fe-S protein 5, 15kDa(NADH-coenzyme Q reductase)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	15 kDa
Uniprot ID	O43920
Dilution Range	WB: 1:1000-1:4000; IHC: 1:20-1200
Expiration Date	12 months from date of receipt.



Western blot analysis of human kidney ti...



Immunohistochemical staining of paraffin...