



Product Datasheet

NF200 antibody (orb395356)

Description	Mouse monoclonal antibody to NF200.
Species/Host	Mouse
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	neurofilament, heavy polypeptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	200 kDa
Uniprot ID	P12036
Dilution Range	WB: 1:1000-1:8000; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of pig brain tissue...



Immunohistochemical staining of paraffin...