



Product Datasheet

NOB1 antibody (orb395364)

Description	Mouse monoclonal antibody to NOB1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	NIN1/RPN12 binding protein 1 homolog(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	47kd
Uniprot ID	Q9ULX3
Dilution Range	WB: 1:1000-1:6000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells us...



Immunohistochemical staining of paraffin...