




Product Datasheet

ATF7 antibody (orb39537)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	ATF7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51757
Uniprot ID	P17544
NCBI	NP_001193612.1 , NP_006847.1 , NP_001123532.1 , NP_001193611.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, FC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of human placenta ...