



Product Datasheet

NUMBL antibody (orb395372)

Description	Mouse monoclonal antibody to NUMBL.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	numb homolog(Drosophila)-like
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	65-70 kDa
Uniprot ID	Q9Y6R0
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of L02 cells using...



Immunohistochemical staining of paraffin...