


Product Datasheet

NY-ESO-1 antibody (orb395374)

Description	Rabbit polyclonal antibody to NY-ESO-1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	cancer/testis antigen 1A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	23 kDa
Uniprot ID	P78358
Dilution Range	WB: 1:200-1:1000; IP: 1:200-1:2000; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells using NY-ESO-1 antibody.