



Product Datasheet

ODZ1 antibody (orb395383)

Description	Rabbit polyclonal antibody to ODZ1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	odz, odd Oz/ten-m homolog 1(Drosophila)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	305kd
Uniprot ID	Q9UKZ4
Dilution Range	WB: 1:1000-1:4000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical staining of paraffin...