



Product Datasheet

p120 Catenin antibody (orb395392)

Description	Mouse monoclonal antibody to p120 Catenin.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	catenin, delta 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	120 kDa
Uniprot ID	O60716
Dilution Range	WB: 1:500-1:2000; IHC: 1:100-1:400
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Immunohistochemical staining of paraffin...