



Product Datasheet

P4HA1 antibody (orb395399)

Description	Mouse monoclonal antibody to P4HA1.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	prolyl 4-hydroxylase, alpha polypeptide I
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	61 kDa
Uniprot ID	P13674
Dilution Range	WB: 1:500-1:2000; IHC: 1:200-1:1000; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of PC-3 cells using P4HA1 antibody.



Immunohistochemical staining of paraffin-embedded tissue sections using P4HA1 antibody.