




Product Datasheet

LOX-1 antibody (orb39542)

Description	Rabbit polyclonal antibody against LOX-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	OLR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30959
Uniprot ID	P78380
NCBI	NP_002534.1 , NP_001166104.1 , NP_001166103.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of HL-60 cell line...



Western blot analysis of 293 cell lysate...